REVIEW

Open Access

Check for

Intravital imaging of orthotopic and ectopic bone

Kunihiko Hashimoto^{1,2}, Takashi Kaito¹, Junichi Kikuta² and Masaru Ishii^{2*}

Abstract

Bone homeostasis is dynamically regulated by a balance between bone resorption by osteoclasts and bone formation by osteoblasts. Visualizing and evaluating the dynamics of bone cells in vivo remain difficult using conventional technologies, including histomorphometry and imaging analysis. Over the past two decades, multiphoton microscopy, which can penetrate thick specimens, has been utilized in the field of biological imaging. Using this innovative technique, the in vivo dynamic motion of bone metabolism-related cells and their interactions has been revealed. In this review, we summarize previous approaches used for bone imaging and provide an overview of current bone tissue imaging methods using multiphoton excitation microscopy.

Keywords: Intravital imaging, Multiphoton microscopy, Regeneration, Bone, Ectopic bone

Background

Bone homeostasis is dynamically regulated internally and externally via interactions among osteoblasts, osteoclasts, osteocytes, other organs, and other cell types. Conventional histological evaluation of bone tissue has been performed on histological sections ex vivo; this evaluation is limited to "static" analysis of cell morphology and gene and protein expression. Bone histomorphometric analysis can quantify the rates of bone formation and resorption during a certain period but cannot visualize real-time "dynamic" behaviors, interactions, and functions among osteoblasts, osteoclasts, osteocytes, and other cell types.

Because conventional laser beams used for microscopic observation cannot penetrate the thick mineralized cortical bone, numerous approaches including micro-computed tomography imaging [1-4], Raman microspectroscopy imaging [5, 6], and magnetic resonance imaging [7] have been utilized for indirectly visualizing the inside of the bone through time-lapsed in vivo imaging. However, it remains difficult to perform real-time analysis of bone dynamics.

Full list of author information is available at the end of the article



The advent of multiphoton microscopy has launched a new era in the field of biological imaging of the bone. The longer wavelength light sources used in multiphoton excitation microscopy enable deeper tissue penetration because less scattering occurs, and there is less laser-induced damage to the tissue; the lower levels of photobleaching of the imaged fluorophores by using near-infrared lasers enables longer observation times compared with those of conventional fluorescent microscopy. Using this innovative technique, the dynamic motion of osteoblasts and osteoclasts and their interactions inside skeletal bone have been clarified. We recently established a novel intravital imaging system for ectopic bone formation as a powerful tool for optimizing bone tissue regeneration. In this review, we provide an overview of current bone tissue imaging including skeletal bone and ectopic bone formation using multiphoton excitation microscopy.

Intravital imaging of skeletal bone with multiphoton excitation microscopy

In vivo imaging can be divided into intravital imaging and tissue explant imaging [8]. In intravital imaging, experimental animals are kept alive under anesthesia, and tissues of interest are exposed for observation. In contrast, in tissue explant imaging, the tissues harvested from experimental animals are maintained under viable

© The Author(s). 2020 **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.org/licenses/by/4.0/.

^{*} Correspondence: mishii@icb.med.osaka-u.ac.jp

²Department of Immunology and Cell Biology, Graduate School of Medicine & Frontier Biosciences, Osaka University, 2-2 Yamada-oka, Suita, Osaka 565-0871, Japan

Table 1 Fluorescent-labeled	mice for intravital mu	ltiphoton imagi	ng in bone	
Target gene/promoter	Fluorescence (color)	Labeled cells	Report (development)	Report intravital imaging)
CX3CR1	EGFP (green)	Macrophage	Jung et al. Mol. Cell. Biol. (2000) [16].	Ishii et al. Nature (2009) [17].
V-ATPase a3 subunit promoter	GFP (green)	Macrophage	Sun-Wada et al. J Cell Sci. (2009) [14].	Kikuta et al. J. Clin. Invest. (2013) [12].
TRAP (exon 2)	tdTomato (red)	Osteoclasts	Kikuta et al. J. Clin. Invest. (2013) [12].	Kikuta et al. J. Clin. Invest. (2013) [12]. Furuya et al. Nat. Commun. (2018) [13].
Rat COL1A1 promoter (2.3 kb)	ECFP (cyan)	Osteoblasts	Furuya et al. Nat. Commun. (2018) [13].	Furuya et al. Nat. Commun. (2018) [13]. Hashimoto et al. Sci Rep. (2020) [15].
Rat COL1A1 promoter (2.3 kb)	GFP (green)	Osteoblasts	Kalajzic et al. J. Bone Miner. Res. (2002) [18].	Kalajzic et al. J. Bone Miner. Res. (2002) [15]. Villa et al. Tissue Eng Part C Methods. (2013) [19]. Huang C, et al. J Bone Miner Res. (2015). [20]
CX3CB1 C-X3-C motif chemokine re-	rentor 1 EGEP enhanced o	ireen fluorescent ni	otein V-4TPase V-type H+-ATPase GEP greep filiore	scent protein 784P tartrate-resistant acid phosphatase COU141 collagen type 1

agen type 5 Ā ÷ v-type escent protein, V-AIPase CX3CR1 C-X3-C motif chemokine receptor 1, EGFP enhanced green fluore alpha 1, ECFP enhanced cyan fluorescent protein

ole 2 Small molecular 1	fluorescent probes for i	ntravital multiphoton imaging	l in bone		
all molecular probe	Fluorescence	Cell function	Report	Report	
			(development)	(intravital imaging)	

Table 2 Small molecular f	luorescent probes for in	ntravital multiphoton imaging in bone		
Small molecular probe	Fluorescence	Cell function	Report (development)	Report (intravital imaging)
BAp-E	Green	Bone resorption (osteoclast)	Kowada et al. J Am Chem Soc. (2011) [21].	Kikuta et al. J. Clin. Invest. (2013) [12].
pHocas-3	Green	Bone resorption (osteoclast)	Maeda et al. Nat. Chem. Biol. (2016) [22].	Furuya et al. Nat. Commun. (2018) [13].

B4p-E boron-dipyrromethene-based (BODIPY-based) H+-sensing fluorescent probe, pHocas-3 pH-activatable fluorescent probe for osteoclast activity sensing-3

Target molecule/protein	Visualization (color)	Target tissue	Report (development)	Report (intravital imaging)
Noncentrosymmetric molecular organization (collagen etc.)	SHG (blue, etc.*)	Bone etc.	Freund et al. Biophys J (1986) [23].	Hashimoto et al. Sci Rep. (2020) [15]. Villa et al. Tissue Eng Part C Methods. (2013) [19].
Ca ²⁺	Alizarin (red)	Bone (mineral)	O'Brien et al. J Biomech. (2002) [26].	Hashimoto et al. Sci Rep. (2020) [15]. Villa et al. Tissue Eng Part C Methods. (2013) [19].

Table 3 Other methods for bone tissue visualization in intravital multiphoton imaging

SHG second harmonic generation

*Color changes based on the wavelength of incidental laser

conditions for observation in culture media. In this review, we discuss only intravital bone imaging with multiphoton excitation microscopy. Mouse calvarial bone, in which the distance from the bone surface to the bone cavity is only 80-120 µm, was used in studies to visualize the bonemarrow cavity. This is because thick cortical bone consists crystalized calcium phosphate, which of prevents visualization at depths greater than 150-200 µm, even with the use of a near-infrared laser [8]. However, in addition to the calvarial bone, there are reports of imaging using multiphoton excitation microscopy on the long bones such as the femur [9] or tibia [10, 11]. Intravital imaging of the marrow of long bones is more invasive because surgical treatment (bone thinning or drilling) of the cortical bone is required. In addition, for the reduction of artifacts (body motion) during the intravital imaging experiments, the imaging box and the anesthetized mouse are maintained at a constant warm temperature using heated air. Heart rate is monitored using an electrocardiogram monitor device, and the concentration of the anesthetic gas is adjusted by using the heart rate as a guide.

Visualization of osteoclasts, osteoblasts, and matrix

To visualize osteoclasts and osteoblasts, genetically modified (knock-in) mice expressing a fluorescent protein in cell-specific genes are employed. To visualize mature osteoclasts, transgenic mice expressing fluorescent proteins under the control of the tartrate-resistant acid phosphatase (TRAP) promoter (TRAP-tdTomato mice) or the ATP-driven proton pump (V-type H+-ATPase a3 subunit) promoter (a3-GFP mice) are used [12–15]. To visualize osteoclast precursors, knock-in mice expressing enhanced green fluorescent protein (EGFP) under the promoter of fractalkine receptor (CX3CR1), which is expressed mainly on monocyte-lineage cells, including osteoclast precursors, was used (CX₃CR1-EGFP mice) [16, 17].

To visualize osteoblasts, transgenic mice expressing GFP or enhanced cyan fluorescent protein (ECFP) under the promoter of type I collagen in osteoblasts (Col 2.3-GFP and Col 2.3-ECFP) is generally used [13, 15, 18–20] (Table 1).

To further visualize the functional status of cells, chemical fluorescent probes are being developed such as pH-activatable probes, by which a low-pH region created by bone-resorbing osteoclasts can be visualized [12, 13, 21, 22] (Table 2).

By using near-infrared lasers for multiphoton excitation, collagen fibers in the bone matrix can be visualized by a nonlinear optical process named as second harmonic generation (SHG), without additional fluorescent labeling [15, 19, 23– 25]. Additionally, administration of alizarin can be performed to visualize calcified bone matrix [15, 19, 26] (Table 3).

Multifluorescent images are acquired via direct detection of fluorescence using four external non-descanned detectors equipped with dichroic and emission filters. The acquired images are subjected to channel unmixing using the NIS-Elements integrated software (Nikon) for autofluorescence and crosstalk reduction [13, 15].

Dynamic visualization of osteoclasts

Osteoclasts, which have a unique hematopoietic origin and move dynamically among bone cells, can be visualized by intravital imaging, which has provided insights into these cells. Osteoclast precursors circulate in the blood vessel and migrate to the bone surface to perform their bone resorptive activity. However, the mechanism underlying the control of osteoclast precursor migration remains unclear. We demonstrated that sphingosine-1-phosphate (S1P), an abundant lipid mediator in the blood, regulates the migration of osteoclast precursors between blood vessels and the bone surface through two S1P receptors (S1PR1 and S1PR2) expressed on osteoclast precursors [17]. The action mechanism of the anti-osteoporosis agent active vitamin D was shown to involve reduced expression of S1PR2 on osteoclast precursors at the bone surface and subsequent migration of osteoclast precursors from the bone surface to the blood vessels [27].

Table 4 Merits and demerits of multiphoton excitation

microscopy
Provide the second sec

Multiphoton excitation micros	сору
Merits	Demerits
1. Time-lapse intravital imaging	1. Visualization of only fluorescently labeled target cells and SHG
2. Imaging without tissue preparation and staining	2. Limitation of penetration depth
3. Visualization of collagen fibers as SHG	3. Difficulty in visualizing inflamed tissues with effusion or tissues with bleeding



Combined visualization and quantification of osteoclast motility and the acidic environment created by osteoclasts using pH-sensing chemical fluorescent probes clarified that mature osteoclasts can be classified into two different types based on their motility and function, that is, the static-bone resorptive (R-type) and movingnon-resorptive (N-type) [12]. Intravital imaging enabled real-time evaluation of the conversion from N-type to Rtype osteoclasts through direct cell-to-cell contact with Th17 cells (a subset of RANKL-expressing CD4+ T cells) and intravenous administration of RANKL to induce the formation of mature osteoclasts.

Dynamic interaction between osteoclasts and osteoblasts

The coupling of bone formation and resorption has been increasingly recognized as a complex, dynamic process. Cell-to-cell contacts between osteoblasts and osteoclasts play an important role in coupling [28], but the spatiotemporal relationship and functional changes caused by the contacts remain unclear. We recently performed intravital imaging of the interaction between osteoblasts and osteoclasts by using double transgenic mice expressing EGFP driven by the type I collagen promoter in osteoblasts and tdTomato under control of the TRAP promoter in osteoclasts [13]. We demonstrated that osteoclasts do not perform bone resorption when they are in contact with osteoblasts and engage in bone resorptive activity when they are not in contact with osteoblasts. These results suggest that osteoblasts control the bone resorption activity of osteoclasts via direct contact.

Intravital imaging of ectopic bone formation

Intravital imaging of skeletal bone by multiphoton excitation microscopy has revealed the in vivo dynamic motion of osteoclasts and osteoblasts, as well as their interactions with each other (Fig. 1). This innovative technique may also be applicable in the field of bone regeneration. To optimize bone regeneration, the process by which new bone is created, it is important to understand the temporospatial appearance and motility of





Fig. 2 Representative intravital two-photon microscopy images of ectopic bone formation area in Col2.3-ECFP mice (visualization of blood vessels). Visualization of osteoblasts, newly formed bone (collagen), and blood vessels after implantation of a collagen sponge containing BMP-2 [15]. Cyan, osteoblasts expressing Col2.3-ECFP; red, blood vessels stained with rhodamine; blue, collagen fibers (second harmonic generation, SHG). Scale bar, 100 μm



osteoclasts and osteoblasts and the bone formation process. Indeed, the investigations on BMSC-mediated calvarial bone defect repair using multiphoton excitation microscopy have been reported [19, 20]. In the field of bone tissue engineering, bone morphogenetic protein (BMP) plays a central role owing to its potent bone induction ability [29]. However, the dosedependent inflammation-related side effects of BMP prevent its widespread use [30, 31]. One possible reason for the difficulty in optimizing BMP-induced bone formation is the limited understanding of this process in vivo. We successfully established a method for intravital imaging of the BMP-2-Induced ectopic bone formation process [15]. The temporospatial dynamic bone formation process, which is initiated by angiogenesis, the recruitment of osteoblasts and osteoclasts, the formation of collagen fibers, and deposition of bone minerals were visualized by intravital imaging of Col2.3-ECFP mice and Col2.3-ECFP/TRAP tdTomato mice (Figs. 2, 3, and 4). Additionally, we detected spindle-shaped osteoblasts and unidirectional collagen fibers, with the osteoblasts and collagen fibers showing an orientational correlation in the early stage of ectopic bone formation (Fig. 5).

Furthermore, it is known that intermittent administration of teriparatide accelerates osteogenic differentiation of mesenchymal progenitor cells instead of adipogenic differentiation [32], reactivates lining cells [33], and suppresses apoptosis in osteoblasts [34]; it also activates bone remodeling, which results in the increased bone



of osteoblasts, osteoclasts, and newly formed bone (collagen) after implantation of a collagen sponge containing BMP-2 [15]. Cyan, osteoblasts expressing Col2.3-ECFP; red, osteoclasts expressing TRAP-tdTomato; blue, collagen fibers (second harmonic generation, SHG). Scale bar, 100 µm



resorptive activity of osteoclasts. Moreover, in this study, intermittent administration of teriparatide was shown to increase BMP-induced bone volume with the increase in the number of osteoblasts, and caused a decrease in the dynamic morphological changes (activated resorptive activity) in osteoclasts.

Intravital imaging can be used to visualize the whole BMP-2-induced bone induction process, contributing to the understanding of ectopic bone formation and providing a foundation for optimizing this process.

Future perspective

Intravital imaging has the greatest advantage with regard to spatiotemporal visualization of living tissues, which cannot be achieved using other methods. However, current imaging techniques have several limitations (Table 4). First, not all objects can be observed in the visual fields of multiphoton microscopy. Because only fluorescent labeling and SHG allow us to see the target cells or organs, we should not misunderstand the lack of a signal as showing an open field. As diverse structures and cellular components are probably present, observations need to be interpreted carefully. Second, the penetration depth in two-photon microscopy is up to 800-1000 µm in soft tissues and 200 µm in hard tissues (e.g., bone). Therefore, there are only limited areas where this imaging can be applied. However, the recently developed three-photon microscopy substantially improves the penetration depth and is expected to overcome these limitations with its resolution [35]. Moreover, because of the wide scattering of light, it is difficult to visualize inflamed tissues with effusion or tissues with bleeding using multiphoton microscopy. To resolve these problems, technical innovations in optical systems and fluorophore are expected in the future.

Conclusions

Intravital imaging has made it possible to visualize the "real-time" dynamics of living cells and biological phenomena caused by cell-to-cell interactions and interventions in various organs. Application of this technique to bone imaging in the past 10 years has provided insights into bone-related processes. Rapid advances in the development of optical instruments and functional chemical fluorescent probes [9] will enable visualization of the dynamics and function in bone.

Abbreviations

TRAP: Tartrate-resistant acid phosphatase; EGFP: Enhanced green fluorescent protein; ECFP: Enhanced cyan fluorescent protein; SHG: Second harmonic generation; S1P: Sphingosine-1-phosphate; BMP: Bone morphogenetic protein; Col 2.3: Promoter of type I collagen in osteoblasts

Acknowledgements

Shigeto Seno and Hideo Matsuda, Department of Bioinformatic Engineering, Graduate School of Information Science & Technology, Osaka University.

Authors' contributions

K.H. and T.K. wrote the initial draft. J.K and M.I. revised the final draft. All authors read and approved the final manuscript.

This work was supported by CREST, Japan Science and Technology Agency; Grants-in-Aid for Scientific Research (S) from the Japan Society for the Promotion of Science (JSPS to M.I.); Grant-in-Aid for Scientific Research (C) from JSPS (to T.K.); Grant-in-Aid for Young Scientists (A) from JSPS (to J.K.); grants from the Uehara Memorial Foundation (to M.I.); from the Kanae Foundation for the Promotion of Medical Sciences (to M.I.); from Mochida Memorial Foundation (to M.I.); and from the Takeda Science Foundation (to M.I.).

Availability of data and materials

Not applicable.

Ethics approval and consent to participate

Not applicable.

Consent for publication

Not applicable.

Competing interests

The authors declare that they have no competing interests.

Author details

¹Department of Orthopedic Surgery, Graduate School of Medicine, Osaka University, Osaka 565-0871, Japan. ²Department of Immunology and Cell Biology, Graduate School of Medicine & Frontier Biosciences, Osaka University, 2-2 Yamada-oka, Suita, Osaka 565-0871, Japan.

Received: 26 May 2020 Accepted: 31 July 2020 Published online: 02 November 2020

References

- Lambers FM, Schulte FA, Kuhn G, Webster DJ, Müller R. Mouse tail vertebrae adapt to cyclic mechanical loading by increasing bone formation rate and decreasing bone resorption rate as shown by time-lapsed in vivo imaging of dynamic bone morphometry. Bone. 2011;49:1340–50.
- Birkhold AI, Razi H, Weinkamer R, Duda GN, Checa S, Willie BM. Monitoring in vivo (re)modeling: A computational approach using 4D microCT data to quantify bone surface movements. Bone. 2015;75:210–21.
- Lienemann PS, Metzger S, Kiveliö AS, Blanc A, Papageorgiou P, Astolfo A, Pinzer BR, Cinelli P, Weber FE, Schibli R, Béhé M, Ehrbar M. Longitudinal *in vivo* evaluation of bone regeneration by combined measurement of multi-pinhole SPECT and micro-CT for tissue engineering. Sci Rep. 2015;5: 10238.
- Lambers FM, Stuker F, Weigt C, Kuhn G, Koch K, Schulte FA, Ripoll J, Rudin M, Müller R. Longitudinal in vivo imaging of bone formation and resorption using fluorescence molecular tomography. Bone. 2013;52:587–95.
- Penela G, Delfossea C, Descampsb M, Leroy G. Composition of bone and apatitic biomaterials as revealed by intravital Raman microspectroscopy. Bone. 2005;36:893–901.
- Sroka-Bartnicka A, Kimber JM, Borkowski L, Pawlowska M, Polkowska I, Kalisz G, Belcarz A, Jozwiak K, Ginalska G, Kazarian SG. The biocompatibility of carbon hydroxyapatite/β-glucan composite for bone tissue engineering studied with Raman and FTIR spectroscopic imaging. Anal Bioanal Chem. 2015;407:7775–85.
- Jing XH, Yang L, Duan XJ, Xie B, Chen W, Li Z, Tan HB. In vivo MR imaging tracking of magnetic iron oxide nanoparticle labeled, engineered, autologous bone marrow mesenchymal stem cells following intra-articular injection. Joint Bone Spine. 2008;75:432–8.
- Ishii M. Intravital imaging technology reveals immune system dynamics in vivo. Allergol Int. 2016;65:225–7.
- Reismann D, Stefanowski J, Günther R, Rakhymzhan A, Matthys R, Nützi R, Zehentmeier S, Schmidt-Bleek K, Petkau G, Chang HD, Naundorf S, Winter Y, Melchers F, Duda G, Hauser AE, Niesner RA. Longitudinal intravital imaging of the femoral bone marrow reveals plasticity within marrow vasculature. Nat Commun. 2017;8:2153.
- Kim S, Lin L, Brown GAJ, Hosaka K, Scott EW. Extended time-lapse in vivo imaging of tibia bone marrow to visualize dynamic hematopoietic stem cell engraftment. Leukemia. 2017;31:1582–92.
- 11. Köhler A, Schmithorst V, Filippi MD, Ryan MA, Daria D, Gunzer M, Geiger H. Altered cellular dynamics and endosteal location of aged early

hematopoietic progenitor cells revealed by time-lapse intravital imaging in long bones. Blood. 2009;114:290–8.

- Kikuta J, Wada Y, Kowada T, Wang Z, Sun-Wada GH, Nishiyama I, Mizukami S, Maiya N, Yasuda H, Kumanogoh A, Kikuchi K, Germain RN, Ishii M. Dynamic visualization of RANKL and Th17-mediated osteoclast function. J Clin Invest. 2013;123:866–73.
- Furuya M, Kikuta J, Fujimori S, Seno S, Maeda H, Shirazaki M, Uenaka M, Mizuno H, Iwamoto Y, Morimoto A, Hashimoto K, Ito T, Isogai Y, Kashii M, Kaito T, Ohba S, Chung UI, Lichtler AC, Kikuchi K, Matsuda H, Yoshikawa H, Ishii M. Direct cell-cell contact between mature osteoblasts and osteoclasts dynamically controls their functions in vivo. Nat Commun. 2018;9:300.
- Sun-Wada GH, Tabata H, Kawamura N, Aoyama M, Wada Y. Direct recruitment of H+-ATPase from lysosomes for phagosomal acidification. J Cell Sci. 2009;122:2504–13.
- Hashimoto K, Kaito T, Furuya M, Seno S, Okuzaki D, Kikuta J, Tsukazaki H, Matsuda H, Yoshikawa H, Ishii M. In vivo dynamic analysis of BMP-2-induced ectopic bone formation. Sci Rep. 2020;10:4751.
- Jung S, Aliberti J, Graemmel P, Sunshine MJ, Kreutzberg GW, Sher A, Littman DR. Analysis of fractalkine receptor CX3CR1 function by targeted deletion and green fluorescent protein reporter gene insertion. Mol Cell Biol. 2000; 20:4106–14.
- Ishii M, Egen JG, Klauschen F, Meier-Schellersheim M, Saeki Y, Vacher J, Proia RL, Germain RN. Sphingosine-1-phosphate mobilizes osteoclast precursors and regulates bone homeostasis. Nature. 2009;458:524–8.
- Kalajzic I, Kalajzic Z, Kaliterna M, Gronowicz G, Clark SH, Lichtler AC, Rowe D. Use of type I collagen green fluorescent protein transgenes to identify subpopulations of cells at different stages of the osteoblast lineage. J Bone Miner Res. 2002;17:15–25.
- Villa MM, Wang L, Huang J, Rowe DW, Wei M. Visualizing osteogenesis in vivo within a cell-scaffold construct for bone tissue engineering using two-photon microscopy. Tissue Eng Part C Methods. 2013;19:839–49.
- Huang C, Ness VP, Yang X, Chen H, Luo J, Brown EB, Zhang X. Spatiotemporal analyses of osteogenesis and angiogenesis via intravital imaging in cranial bone defect repair. J Bone Miner Res. 2015; 30:1217–30.
- Kowada T, Kikuta J, Kubo A, Ishii M, Maeda H, Mizukami S, Kikuchi K. In vivo fluorescence imaging of bone-resorbing osteoclasts. J Am Chem Soc. 2011; 133:17772–6.
- Maeda H, Kowada T, Kikuta J, Furuya M, Shirazaki M, Mizukami S, Ishii M, Kikuchi K. Real-time intravital imaging of pH variation associated with osteoclast activity. Nat Chem Biol. 2016;12:579–85.
- Freund I, Deutsch M, Sprecher A. Connective tissue polarity. Optical secondharmonic microscopy, crossed-beam summation, and small-angle scattering in rat-tail tendon. Biophys J. 1986;50:693–712.
- 24. Helmchen F, Denk W. Deep tissue two-photon microscopy. Nat Methods. 2005. 2:932–40.
- Campagnola PJ, Loew LM. Second-harmonic imaging microscopy for visualizing biomolecular arrays in cells, tissues and organisms. Nat Biotechnol. 2003;21:1356–60.
- O'Brien FJ, Taylor D, Lee TC. An improved labelling technique for monitoring microcrack growth in compact bone. J Biomech. 2002;35:523–6.
- Kikuta J, Kawamura S, Okiji F, Shirazaki M, Sakai S, Saito H, Ishii M. Sphingosine-1-phosphate-mediated precursor monocyte migration is a critical point of control in antibone-resorptive action of active vitamin D. Proc Natl Acad Sci U S A. 2013;110:7009–13.
- Sims NA, Martin TJ. Coupling the activities of bone formation and resorption: a multitude of signals within the basic multicellular unit. Bonekey Rep. 2014;3:481.
- Wang RN, Green J, Wang Z, Deng Y, Qiao M, Peabody M, Zhang Q, Ye J, Yan Z, Denduluri S, Idowu O, Li M, Shen C, Hu A, Haydon RC, Kang R, Mok J, Lee MJ, Luu HL, Shi LL. Bone morphogenetic protein (BMP) signaling in development and human diseases. Genes Dis. 2014;1: 87–105.
- Guzman JZ, Merrill RK, Kim JS, Overley SC, Dowdell JE, Somani S, Hecht AC, Cho SK, Qureshi SA. Bone morphogenetic protein use in spine surgery in the United States: how have we responded to the warnings? Spine J. 2017; 17:1247–54.
- Tannoury CA, An HS. Complications with the use of bone morphogenetic protein 2 (BMP-2) in spine surgery. Spine J. 2014;14:552–9.
- Balani DH, Ono N, Kronenberg HM. Parathyroid hormone regulates fates of murine osteoblast precursors in vivo. J Clin Invest. 2017;127:3327–38.

- Kim SW, et al. Intermittent parathyroid hormone administration converts quiescent lining cells to active osteoblasts. J Bone Miner Res. 2012;27: 2075–84.
- Jilka RL, Weinstein RS, Bellido T, Roberson P, Parfitt AM, Manolagas SC, et al. Increased bone formation by prevention of osteoblast apoptosis with parathyroid hormone. J. Clin. Invest. 1999;104:439–46.
- 35. Kim JM, Bixel MG. Intravital multiphoton imaging of the bone and bone marrow environment. Cytometry A. 2020;97:496–503.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

